




Pusch®

There is no substitute for quality

KR Version



Version KR belongs to our series of instruments with passive safety. Rounded edges are a characteristic feature of these instruments. With these, nails can be smoothed and callus removed gently and with maximum safety.

New!

New!

854KR

841KR

840KR






Our series of instruments with passive safety is specially designed for safe everyday use by pedicurists and podiatrists. The primary aim, particularly when treating diabetic patients, is to avoid injury at all times. Our series of instruments with passive safety assists the user in this.

The upper edge of the working part is rounded on all instruments of version KR. This ensures gentle treatment of high-risk patients, whether when smoothing nails or removing callus. In addition to the time-proven cylindrical basic forms, we are now offering the smaller cylindrical form 841KR for more delicate treatment. The conical version 854KR is also new: it is available in two sizes so that you can always choose the right instrument for the application and size of nail concerned. Recommended speed: 5,000 – 10,000 rpm.








As a matter of principle, we recommend that the speed be reduced when working with the dry technique, as well as when treating high-risk patients.










Rounded edge



Medium grit

						
BUSCH	841KR	New!	840KR	New!	New!	854KR
Ø ISO	035		055			040

Fine grit

						
BUSCH	8841KR	New!	8840KR	New!	New!	8854KR
Ø ISO	035		055			040

 = Packaging unit



BUSCH & CO. GmbH & Co. KG

Postfach 1152 • 51751 Engelskirchen • Unterkaltenbach 17-27 • 51766 Engelskirchen • GERMANY
 Telefon +49 (0) 2263 860 • Telefax +49 (0) 2263 20741 • mail@busch.eu • www.busch.eu